

RECEIVED

Access DB# 7595E

## SEARCH REQUEST FORM

SEP 17 2002

Scientific and Technical Information Center

STIC

Requester's Full Name: Lynda Guo Examiner #: 79756 Date: 09/17/02  
Art Unit: 1627 Phone Number 30-605-1200 Serial Number: 09/806,983  
Mail Box and Bldg/Room Location: CM1-3B19 Results Format Preferred (circle): PAPER DISK E-MAIL  
office - 3D08

If more than one search is submitted, please prioritize searches in order of need.

\*\*\*\*\*

Please provide a detailed statement of the search topic, and describe as specifically as possible the subject matter to be searched. Include the elected species or structures, keywords, synonyms, acronyms, and registry numbers, and combine with the concept or utility of the invention. Define any terms that may have a special meaning. Give examples or relevant citations, authors, etc, if known. Please attach a copy of the cover sheet, pertinent claims, and abstract.

Title of Invention: Method for determining alkaline phosphatase

Inventors (please provide full names): \_\_\_\_\_

Earliest Priority Filing Date: \_\_\_\_\_

\*For Sequence Searches Only\* Please include all pertinent information (parent, child, divisional, or issued patent numbers) along with the appropriate serial number.

Point of Contact:  
Thomas G. Larson, Ph.D.  
703-308-7309  
CM1, Rm. 6 B 01

## STAFF USE ONLY

Searcher: Point of Contact:  
Thomas G. Larson, Ph.D.  
Searcher Phone #: 703-308-7309  
Searcher Location: CM1, Rm. 6 B 01

Date Searcher Picked Up: 9/20  
Date Completed: 9/27/02  
Searcher Prep & Review Time: 30

Clerical Prep Time: \_\_\_\_\_

Online Time: 289

## Type of Search

NA Sequence (#) \_\_\_\_\_ STN 4391  
AA Sequence (#) \_\_\_\_\_ Dialog \_\_\_\_\_  
Structure (#) \_\_\_\_\_ Questel/Orbit \_\_\_\_\_  
Bibliographic X Dr. Link \_\_\_\_\_  
Litigation \_\_\_\_\_ Lexis/Nexis \_\_\_\_\_  
Fulltext \_\_\_\_\_ Sequence Systems \_\_\_\_\_  
Patent Family \_\_\_\_\_ WWW/Internet \_\_\_\_\_  
Other \_\_\_\_\_ Other (specify) \_\_\_\_\_

## Vendors and cost where applicable

=> d que 138

L34 139 SEA FILE=BIOSIS ABB=ON PLU=ON 9001-78-9# AND 330-13-2#  
L37 25086 SEA FILE=BIOSIS ABB=ON PLU=ON BLOOD/CT  
L38 0 SEA FILE=BIOSIS ABB=ON PLU=ON L37 AND L34

=> d que 140

L34 139 SEA FILE=BIOSIS ABB=ON PLU=ON 9001-78-9# AND 330-13-2#  
L39 10 SEA FILE=BIOSIS ABB=ON PLU=ON NITROPHENYL PHOSPHATE/TI AND  
ALKALINE PHOSPHATASE/TI  
L40 6 SEA FILE=BIOSIS ABB=ON PLU=ON L34 AND L39

=> file biotechno

FILE 'BIOTECHNO' ENTERED AT 16:58:20 ON 26 SEP 2002

COPYRIGHT (C) 2002 Elsevier Science B.V., Amsterdam. All rights reserved.

FILE LAST UPDATED: 24 SEP 2002 <20020924/UP>

FILE COVERS 1980 TO DATE.

>>> SIMULTANEOUS LEFT AND RIGHT TRUNCATION AVAILABLE IN  
/CT AND BASIC INDEX <<<

=> d que 147

L41 36 SEA FILE=BIOTECHNO ABB=ON PLU=ON 9001-78-9# AND 330-13-2#  
L42 304217 SEA FILE=BIOTECHNO ABB=ON PLU=ON BLOOD OR PLASMA OR SERUM  
L43 8 SEA FILE=BIOTECHNO ABB=ON PLU=ON L41 AND L42  
L46 1549 SEA FILE=BIOTECHNO ABB=ON PLU=ON ALKALINE PHOSPHATASE/TI  
L47 3 SEA FILE=BIOTECHNO ABB=ON PLU=ON L46 AND L43

=> d que 149

L41 36 SEA FILE=BIOTECHNO ABB=ON PLU=ON 9001-78-9# AND 330-13-2#  
L42 304217 SEA FILE=BIOTECHNO ABB=ON PLU=ON BLOOD OR PLASMA OR SERUM  
L43 8 SEA FILE=BIOTECHNO ABB=ON PLU=ON L41 AND L42  
L48 6848 SEA FILE=BIOTECHNO ABB=ON PLU=ON SPECTROPHOTOM?  
L49 2 SEA FILE=BIOTECHNO ABB=ON PLU=ON L43 AND L48

=> d que 152

L46 1549 SEA FILE=BIOTECHNO ABB=ON PLU=ON ALKALINE PHOSPHATASE/TI  
L52 2 SEA FILE=BIOTECHNO ABB=ON PLU=ON L46 AND NITROPHENYL  
PHOSPHATE/TI

=> s 147 or 149 or 152

L55 6 L47 OR L49 OR L52

=> file medline

FILE 'MEDLINE' ENTERED AT 13:26:23 ON 27 SEP 2002

FILE LAST UPDATED: 26 SEP 2002 (20020926/UP). FILE COVERS 1958 TO DATE.

On June 9, 2002, MEDLINE was reloaded. See HELP RLOAD for details.

MEDLINE thesauri in the /CN, /CT, and /MN fields incorporate the  
MeSH 2002 vocabulary. Enter HELP THESAURUS for details.

THIS FILE CONTAINS CAS REGISTRY NUMBERS FOR EASY AND ACCURATE

## SUBSTANCE IDENTIFICATION.

=&gt; d que 158

L56	1579	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	9001-78-9#
L57	288	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	330-13-2#
L58	0	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	L56 AND L57

} Medline is not  
very good about  
including CAS  
registry numbers  
in their indexing

=&gt; d que 167

L56	1579	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	9001-78-9#
L57	288	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	330-13-2#
L59	34205	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	ALKALINE PHOSPHATASE+PFT/CT
L60	35784	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	L56 OR L59
L61	11386	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	NITROPHENOLS+PFT/CT
L62	66050	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	ORGANOPHOSPHORUS COMPOUNDS+PFT/CT
L63	623	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	L61 AND L62
L64	1522	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	L63 OR L57 OR ((4 OR P) (W) NITROPHENYL (W) PHOSPHATE)
L65	251	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	L60 AND L64
L66	66607	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	BLOOD CHEMICAL ANALYSIS+NT, PFT/CT
L67	2	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	L65 AND L66

=&gt; d que 178

L57	288	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	330-13-2#
L59	34205	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	ALKALINE PHOSPHATASE+PFT/CT
L61	11386	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	NITROPHENOLS+PFT/CT
L62	66050	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	ORGANOPHOSPHORUS COMPOUNDS+PFT/CT
L63	623	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	L61 AND L62
L64	1522	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	L63 OR L57 OR ((4 OR P) (W) NITROPHENYL (W) PHOSPHATE)
L70	94484	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	SPECTROPHOTOMETRY+NT, PFT/CT
L75	3615	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	L59/MAJ (L) BL/CT
L77	35	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	L64 AND L75
L78	5	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	L77 AND L70

=&gt; d que 180

L56	1579	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	9001-78-9#
L57	288	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	330-13-2#
L59	34205	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	ALKALINE PHOSPHATASE+PFT/CT
L60	35784	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	L56 OR L59
L61	11386	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	NITROPHENOLS+PFT/CT
L62	66050	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	ORGANOPHOSPHORUS COMPOUNDS+PFT/CT
L63	623	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	L61 AND L62
L64	1522	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	L63 OR L57 OR ((4 OR P) (W) NITROPHENYL (W) PHOSPHATE)
L65	251	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	L60 AND L64
L79	16969	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	HEMOLYSIS+PFT/CT
L80	2	SEA	FILE=MEDLINE	ABB=ON	PLU=ON	L65 AND L79

=&gt; s 185 or 167 or 178 or 180

L81 8 S L58 OR L67 OR L78 OR L80

=&gt; file wpids

FILE 'WPIDS' ENTERED AT 13:27:55 ON 27 SEP 2002  
COPYRIGHT (C) 2002 THOMSON DERWENT

FILE LAST UPDATED: 24 SEP 2002 <20020924/UP>  
MOST RECENT DERWENT UPDATE 200261 <200261/DW>  
DERWENT WORLD PATENTS INDEX SUBSCRIBER FILE, COVERS 1963 TO DATE

*No controlled term  
indexing or CAS  
Registry No. indexing  
available in  
WPIDS - search  
by free text.*

>>> The BATCH option for structure searches has been  
enabled in WPINDEX/WPIDS and WPIX >>>

>>> PATENT IMAGES AVAILABLE FOR PRINT AND DISPLAY >>>

>>> FOR DETAILS OF THE PATENTS COVERED IN CURRENT UPDATES,  
SEE <http://www.derwent.com/dwpi/updates/dwpicov/index.html> <<<

>>> FOR A COPY OF THE DERWENT WORLD PATENTS INDEX STN USER GUIDE,  
PLEASE VISIT:  
[http://www.stn-international.de/training\\_center/patents/stn\\_guide.pdf](http://www.stn-international.de/training_center/patents/stn_guide.pdf) <<<

>>> FOR INFORMATION ON ALL DERWENT WORLD PATENTS INDEX USER  
GUIDES, PLEASE VISIT:  
[http://www.derwent.com/userguides/dwpi\\_guide.html](http://www.derwent.com/userguides/dwpi_guide.html) <<<

=> d que 184

L82 77 SEA FILE=WPIDS ABB=ON PLU=ON ALKALINE PHOSPHATASE (5A)  
(SERUM OR BLOOD OR PLASMA)  
L83 94 SEA FILE=WPIDS ABB=ON PLU=ON ((4 OR P) (W) NITROPHENYL (W)  
PHOSPHATE)  
L84 4 SEA FILE=WPIDS ABB=ON PLU=ON L82 AND L83

=> d que 187

L83 94 SEA FILE=WPIDS ABB=ON PLU=ON ((4 OR P) (W) NITROPHENYL (W)  
PHOSPHATE)  
L85 6283 SEA FILE=WPIDS ABB=ON PLU=ON HEMOGLOB? OR HAEMOGLOB?  
L86 33 SEA FILE=WPIDS ABB=ON PLU=ON L85 AND ALKALINE PHOSPHATASE  
L87 2 SEA FILE=WPIDS ABB=ON PLU=ON L86 AND L83

=> d que 190

L88 20621 SEA FILE=WPIDS ABB=ON PLU=ON SPECTROSCOP? OR SPECTROPHOTOM?  
L89 25 SEA FILE=WPIDS ABB=ON PLU=ON L88 AND ALKALINE PHOSPHATASE  
L90 1 SEA FILE=WPIDS ABB=ON PLU=ON L89 AND NITROPHENYL PHOSPHATE

=> d que 193

L91 1326 SEA FILE=WPIDS ABB=ON PLU=ON HEMOLYS? OR HAEMOLYS?  
L92 3 SEA FILE=WPIDS ABB=ON PLU=ON L91 AND ALKALINE PHOSPHATASE  
L93 0 SEA FILE=WPIDS ABB=ON PLU=ON L92 AND NITROPHENYL PHOSPHATE

=> s 184 or 187 or 190 or 193

L94 6 S L84 OR L87 OR L90 OR L93

=> dup rem 181 153 140 155 154 194

FILE 'MEDLINE' ENTERED AT 13:46:54 ON 27 SEP 2002

*remove duplicate documents*

FILE 'CAPLUS' ENTERED AT 13:46:54 ON 27 SEP 2002

USE IS SUBJECT TO THE TERMS OF YOUR STN CUSTOMER AGREEMENT.

PLEASE SEE "HELP USAGETERMS" FOR DETAILS.  
COPYRIGHT (C) 2002 AMERICAN CHEMICAL SOCIETY (ACS)

FILE 'BIOSIS' ENTERED AT 13:46:54 ON 27 SEP 2002  
COPYRIGHT (C) 2002 BIOLOGICAL ABSTRACTS INC.(R)

FILE 'BIOTECHNO' ENTERED AT 13:46:54 ON 27 SEP 2002  
COPYRIGHT (C) 2002 Elsevier Science B.V., Amsterdam. All rights reserved.

FILE 'EMBASE' ENTERED AT 13:46:54 ON 27 SEP 2002  
COPYRIGHT (C) 2002 Elsevier Science B.V. All rights reserved.

FILE 'WPIDS' ENTERED AT 13:46:54 ON 27 SEP 2002  
COPYRIGHT (C) 2002 THOMSON DERWENT  
PROCESSING COMPLETED FOR L81  
PROCESSING COMPLETED FOR L53  
PROCESSING COMPLETED FOR L40  
PROCESSING COMPLETED FOR L55  
PROCESSING COMPLETED FOR L54  
PROCESSING COMPLETED FOR L94  
L96 32 DUP REM L81 L53 L40 L55 L54 L94 (5 DUPLICATES REMOVED)

=> D IBIB ABS 1-32

L96 ANSWER 1 OF 32 CAPLUS COPYRIGHT 2002 ACS  
ACCESSION NUMBER: 2000:238439 CAPLUS  
DOCUMENT NUMBER: 132:247995  
TITLE: Spectrophotometric method for the determination of  
**alkaline phosphatase** in blood serum  
that eliminates the interference of hemoglobin  
INVENTOR(S): Weisheit, Ralph; Treiber, Wolfgang  
PATENT ASSIGNEE(S): Roche Diagnostics G.m.b.H., Germany  
SOURCE: Ger. Offen., 6 pp.  
CODEN: GWXXBX  
DOCUMENT TYPE: Patent  
LANGUAGE: German  
FAMILY ACC. NUM. COUNT: 1  
PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
DE 19846301	A1	20000413	DE 1998-19846301	19981008
WO 2000022162	A1	20000420	WO 1999-EP7394	19991005
W: AE, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, CA, CH, CN, CU, CZ, DE, DK, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, UA, UG, US, UZ, VN, YU, ZA, ZW, AM, AZ, BY, KG, KZ, MD, RU, TJ, TM				
RW: GH, GM, KE, LS, MW, SD, SL, SZ, TZ, UG, ZW, AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG				
AU 9964682	A1	20000501	AU 1999-64682	19991005
EP 1119641	A1	20010801	EP 1999-952500	19991005
R: AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE, SI, LT, LV, FI, RO				
JP 2002527076	T2	20020827	JP 2000-576052	19991005
PRIORITY APPLN. INFO.: DE 1998-19846301 A 19981008 WO 1999-EP7394 W 19991005				

AB The invention concerns the spectrophotometric detn. of alk. phosphatase in blood serum or plasma at two wavelengths thus eliminating the interference of Hb or blood substitutes. Blood substitutes are modified Hb, bovine Hb, recombinant Hb. Hb content can be up to 6500 mg/dL. The selected wavelengths are 450 nm in combination with 480 nm, 546 nm, or 575 nm. Thus 4-nitrophenylphosphate was used as substrate; absorption was measured within 4 min.

L96 ANSWER 2 OF 32 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 2000:238438 CAPLUS  
 DOCUMENT NUMBER: 132:247994  
 TITLE: Spectrophotometric method for the determination of **alkaline phosphatase** in plasma or serum using the rate-blank method to eliminate hemoglobin interference  
 INVENTOR(S): Weisheit, Ralph; Treiber, Wolfgang  
 PATENT ASSIGNEE(S): Roche Diagnostics G.m.b.H., Germany  
 SOURCE: Ger. Offen., 6 pp.  
 CODEN: GWXXBX  
 DOCUMENT TYPE: Patent  
 LANGUAGE: German  
 FAMILY ACC. NUM. COUNT: 1  
 PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
DE 19846300	A1	20000413	DE 1998-19846300	19981008
WO 2000022161	A1	20000420	WO 1999-EP7366	19991005
W: AE, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, CA, CH, CN, CU, CZ, DE, DK, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, UA, UG, US, UZ, VN, YU, ZA, ZW, AM, AZ, BY, KG, KZ, MD, RU, TJ, TM				
RW: GH, GM, KE, LS, MW, SD, SL, SZ, TZ, UG, ZW, AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG				
AU 9961997	A1	20000501	AU 1999-61997	19991005
EP 1119640	A1	20010801	EP 1999-948931	19991005
R: AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT, IE, SI, LT, LV, FI, RO				
JP 2002527075	T2	20020827	JP 2000-576051	19991005
PRIORITY APPLN. INFO.: DE 1998-19846300 A 19981008				
WO 1999-EP7366 W 19991005				

AB The invention concerns the spectrophotometric detn. of alk. phosphatase in blood plasma or serum at 450 nm and 660 nm in combination with the rate-blank method to eliminate the interference of Hb or blood substitutes. To measure the rate of the blank, absorption is measured without adding the 4-nitrophenyl phosphate substrate to the samples.

L96 ANSWER 3 OF 32 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 2000:662602 CAPLUS  
 DOCUMENT NUMBER: 134:14572  
 TITLE: Mechanism of interference of a polymerized hemoglobin blood substitute in an alkaline phosphatase method  
 AUTHOR(S): Chance, Jeffrey J.; Norris, Edward J.; Kroll, Martin H.  
 CORPORATE SOURCE: Division of Clinical Chemistry, Department of Pathology, Johns Hopkins University School of Medicine, Baltimore, MD, 21287-7065, USA

SOURCE: *of interest*  
*A for ref. list*  
 Clinical Chemistry (Washington, D. C.) (2000), 46(9), 1331-1337  
 CODEN: CLCHAU; ISSN: 0009-9147  
 PUBLISHER: American Association for Clinical Chemistry  
 DOCUMENT TYPE: Journal  
 LANGUAGE: English

AB Background: Hb-based oxygen carriers can cause profound interferences in many anal. procedures. We detd. the mechanism of interference in the assay of alk. phosphatase activity and identified approaches that might be used to correct for this interference. Methods: Interference of a polymd. Hb blood substitute with the assay of alk. phosphatase was examd. with a Hitachi 917 analyzer and UV-visible spectrophotometry. Results: Hb-based oxygen carrier solns. had substantial absorbance at 415 nm, the wavelength of anal. used to measure the formation of 4-nitrophenol. In addn. to offsetting the initial absorbance at the anal. wavelength, polymd. Hb gave rise to a strong neg. interference plot because of alkali denaturation of the substitute. The same interference mechanism was also obsd. for native Hb (hemolyzate), indicating that the interference was not derived from the polymn. process. The interference can be cor. by implementing a rate-correction procedure, or the interference can be avoided by measurement at 450 nm. Conclusions: The interference of polymd. Hb in the alk. phosphatase assay is a result of an absorbance offset caused by alkali denaturation of Hb. The interference can be cor. or avoided by modifying the calcn. or the anal. wavelength. The correction strategy may also be applicable to improving the hemolysis index for this method.

REFERENCE COUNT: 11 THERE ARE 11 CITED REFERENCES AVAILABLE FOR THIS RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

L96 ANSWER 4 OF 32 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 1999:487468 CAPLUS

DOCUMENT NUMBER: 131:127388

TITLE: Detection system using liposomes and signal modification

INVENTOR(S): Nicklin, Stephen; Clarke, David John; Lloyd, Christopher James; Aojula, Harmesh Singh; Tsilosani, Marina; Wilson, Michael Thomas

PATENT ASSIGNEE(S): The Secretary of State for Defence, UK

SOURCE: PCT Int. Appl., 111 pp.

CODEN: PIXXD2

DOCUMENT TYPE: Patent

LANGUAGE: English

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
WO 9938009	A1	19990729	WO 1999-GB208	19990121
W:	AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, CA, CH, CN, CU, CZ, DE, DK, EE, ES, FI, GB, GD, GE, GH, GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR, LS, LT, LU, LV, MD, MG, MK, MN, MW, MX, NO, NZ, PL, PT, RO, RU, SD, SE, SG, SI, SK, SL, TJ, TM, TR, TT, UA, UG, US, UZ, VN, YU, ZW, AM, AZ, BY, KG, KZ, MD, RU, TJ, TM			
RW:	GH, GM, KE, LS, MW, SD, SZ, UG, ZW, AT, BE, CH, CY, DE, DK, ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, BF, BJ, CF, CG, CI, CM, GA, GN, GW, ML, MR, NE, SN, TD, TG			
ZA 9900325	A	19990719	ZA 1999-325	19990118
CA 2318170	AA	19990729	CA 1999-2318170	19990121
AU 9921770	A1	19990809	AU 1999-21770	19990121
AU 749955	B2	20020704		

EP 1049932 A1 20001108 EP 1999-901770 19990121  
 R: AT, BE, CH, DE, DK, ES, FR, GB, IT, LI, LU, NL, SE, MC, PT, IE, FI  
 NO 200003709 A 20000921 NO 2000-3709 20000719  
 PRIORITY APPLN. INFO.: GB 1998-1120 A 19980121  
 WO 1999-GB208 W 19990121

AB A process for detecting an analyte comprises (a) contacting a sample suspected of contg. said analyte with a containment means comprising a barrier which separates signal generating reagents from said sample, in the presence of an element which interacts specifically with said analyte, under conditions whereby interaction between the analyte and the said element results in activation of the signal generating reagents within the containment means on the side of the barrier opposite to the sample, and (b) detecting any signal generated and retained within the containment means from the sample side of the barrier. The process of the invention provides for sensitive detection of very small nos. of analyte materials using measurement techniques which include counting methods such as flow cytometry. TNT was detected using Tris-HCl pH 7.4, TNP-conjugated melittin as pore-forming reagent, liposomes contg. alk. phosphatase, ELF-97 substrate, and monoclonal antibodies to TNT. Fluorescent liposomes were counted.

REFERENCE COUNT: 7 THERE ARE 7 CITED REFERENCES AVAILABLE FOR THIS RECORD. ALL CITATIONS AVAILABLE IN THE RE FORMAT

L96 ANSWER 5 OF 32 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 1999:431645 CAPLUS  
 DOCUMENT NUMBER: 131:99250  
 TITLE: A stabilized reagent for measuring **alkaline phosphatase** activity.  
 INVENTOR(S): Kumasawa, Ichiro  
 PATENT ASSIGNEE(S): Nissho Corp., Japan  
 SOURCE: Jpn. Kokai Tokkyo Koho, 4 pp.  
 CODEN: JKXXAF  
 DOCUMENT TYPE: Patent  
 LANGUAGE: Japanese  
 FAMILY ACC. NUM. COUNT: 1  
 PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
JP 11187897	A2	19990713	JP 1997-357330	19971225

AB A method is described for improving the storage stability of p-nitrophenylphosphate soln. used as a substrate for measuring alk. phosphatase (ALP) activity in a biol. sample (e.g., blood, urine). In order to suppress the hydrolytic decompn. of p-nitrophenylphosphate upon storage, a substance possessing at least one carboxylic acid group or its salt is added to the soln. of p-nitrophenylphosphoric acid or its salt. The substrate reagent prepd. by this method for measuring ALP provides assay results excellent in reproducibility without giving any inhibitory effect on ALP reaction.

L96 ANSWER 6 OF 32 EMBASE COPYRIGHT 2002 ELSEVIER SCI. B.V.

ACCESSION NUMBER: 1999125970 EMBASE  
 TITLE: Isotachophoretic determination of phosphate splitting from Amifostine and p-nitrophenyl phosphate in serum and neuroblastoma cells.  
 AUTHOR: Renner S.; Klingebiel T.; Niethammer D.; Bruchelt G.; Meissner T.; Eisenbeiss F.  
 CORPORATE SOURCE: S. Renner, Children's Hospital, University of Tübingen, Rumelinstrasse 23, D-72070 Tübingen, Germany.  
 chrrenn@aol.com



SOURCE: Journal of Chromatography A, (1999) 838/1-2 (251-257).  
Refs: 9  
ISSN: 0021-9673 CODEN: JCRAEY  
PUBLISHER IDENT.: S 0021-9673(99)00084-9  
COUNTRY: Netherlands  
DOCUMENT TYPE: Journal; Conference Article  
FILE SEGMENT: 016 Cancer  
029 Clinical Biochemistry  
037 Drug Literature Index  
LANGUAGE: English  
SUMMARY LANGUAGE: English  
AB Amifostine [WR-2721; H<sub>2</sub>N-(CH<sub>2</sub>)<sub>3</sub>-NH-(CH<sub>2</sub>)<sub>2</sub>-S-PO<sub>3</sub>H<sub>2</sub>] is used as a protecting agent in the chemotherapy of neuroblastoma. It is supposed that Amifostine will be transformed into its active form, the free thiol (WR-1065), easier by normal cells than by tumour cells. Analytical capillary isotachophoresis was used to determine the dephosphorylation of Amifostine in serum and on neuroblastoma cells and peripheral blood cells. Furthermore, the biological effects of Amifostine and its free thiol, on cell proliferation of neuroblastoma cells were measured in combination with Carboplatin. It was found that neuroblastoma cells did not split phosphate less efficiently than normal peripheral blood cells. Furthermore, neither Amifostine (as expected) nor the free thiol (not expected according to the theory) were able to inhibit the effects of Carboplatin. Therefore, the current hypothesis concerning the mode of action of Amifostine must be questioned. Copyright (C) 1999 Elsevier Science B.V.

L96 ANSWER 7 OF 32 BIOTECHNO COPYRIGHT 2002 Elsevier Science B.V.DUPLICATE  
ACCESSION NUMBER: 1999:29238108 BIOTECHNO  
TITLE: PCR-ELISA for diagnosis of mucocutaneous leishmaniasis  
AUTHOR: Pinero J.; Martinez E.; Pacheco R.; Aragon Z.; De Armas F.; Del Castillo A.; Valladares B.  
CORPORATE SOURCE: J. Pinero, Departamento Parasitologia, Facultad de Farmacia, Universidad de la Laguna, Avda. Astrofisico Fco. Sanchez s.n., Tenerife, Islas Canarias, Spain.  
E-mail: jpinero@ull.es  
SOURCE: Acta Tropica, (1999), 73/1 (21-29), 16 reference(s)  
CODEN: ACTRAQ ISSN: 0001-706X  
PUBLISHER ITEM IDENT.: S0001706X99000157  
DOCUMENT TYPE: Journal; Article  
COUNTRY: Netherlands  
LANGUAGE: English  
SUMMARY LANGUAGE: English

AN 1999:29238108 BIOTECHNO  
AB In this work we demonstrate that the PCR-ELISA technique is sufficiently sensitive and specific for use as a diagnostic test in cases of mucocutaneous leishmaniasis. DNA was extracted from cultures of *Leishmania braziliensis*, *Leishmania infantum*, *Leishmania tropica*, *Leishmania mexicana*, *Trypanosoma cruzi*, and **blood** samples from individuals who presented a clinical diagnosis of leishmaniasis as well as from healthy individuals. The DNA was PCR amplified and the product obtained was hybridised with a biotin-labelled probe, the sequence of which was designed in our laboratory. The result of the hybridisation was visualised by means of an ELISA technique using antiluorescein antibody labelled with alkaline phosphatase and p-nitrophenylphosphate (pNFF) as chromogen. The optical density of the products of the pNFF hydrolysis was quantified in a **spectrophotometer** at a wavelength of 405 nm. Using this technique the percentage of detection was 83.3% in **blood** samples from patients clinically diagnosed as having mucocutaneous leishmaniasis. No false positive results were obtained.

Copyright (C) 1999 Elsevier Science B.V.

L96 ANSWER 8 OF 32 WPIDS (C) 2002 THOMSON DERWENT  
 ACCESSION NUMBER: 1998-460115 [40] WPIDS  
 DOC. NO. CPI: C1998-139114  
 TITLE: Dry analysis component for **alkaline phosphatase** assay in **blood serum** - includes fine particles of amine salt of **p-nitrophenyl phosphate** or non-crystalline fine particle of amino alcohol salt.  
 DERWENT CLASS: B04 D16  
 PATENT ASSIGNEE(S): (FUJF) FUJI PHOTO FILM CO LTD  
 COUNTRY COUNT: 1  
 PATENT INFORMATION:

PATENT NO	KIND	DATE	WEEK	LA	PG
JP 10191997	A	19980728	(199840)*		6

## APPLICATION DETAILS:

PATENT NO	KIND	APPLICATION	DATE
JP 10191997	A	JP 1997-17919	19970116

PRIORITY APPLN. INFO: JP 1997-17919 19970116  
 AN 1998-460115 [40] WPIDS  
 AB JP 10191997 A UPAB: 19981008  
 Dry analysis component for **alkaline phosphatase** assay in **blood serum** contains fine particles of amine salt of **p-nitrophenyl phosphate** or non-crystalline fine particles of amino alcohol salt in a dispersion.  
 ADVANTAGE - The component gives good reproducibility of analysis value by quantitative analysis of alkaline phosphate. Analysis results with high accuracy are obtained.  
 Dwg.0/0

L96 ANSWER 9 OF 32 BIOSIS COPYRIGHT 2002 BIOLOGICAL ABSTRACTS INC.  
 ACCESSION NUMBER: 1997:334949 BIOSIS  
 DOCUMENT NUMBER: PREV199799634152  
 TITLE: Evaluation of various liquid **p-nitrophenyl phosphate** (pNPP) substrates of **alkaline phosphatase**.  
 AUTHOR(S): Cook, N. (1); Macinnes, B. (1); Rajadhyaksha, M. (1); Kumar, V.  
 CORPORATE SOURCE: (1) IMMCO Diagnostics Inc., Buffalo, NY USA  
 SOURCE: Clinical Chemistry, (1997) Vol. 43, No. 6 PART 2, pp. S249. Meeting Info.: 49th Annual Meeting of the American Association for Clinical Chemistry Atlanta, Georgia, USA July 20-24, 1997  
 ISSN: 0009-9147.  
 DOCUMENT TYPE: Conference; Abstract; Conference  
 LANGUAGE: English

L96 ANSWER 10 OF 32 BIOTECHNO COPYRIGHT 2002 Elsevier Science B.V.DUPLICATE  
 ACCESSION NUMBER: 1995:25337878 BIOTECHNO  
 TITLE: Monoclonal antibody assay for measuring bone-specific **alkaline phosphatase** activity in **serum**

AUTHOR: Gomez Jr. B.; Ardakani S.; Ju J.; Jenkins D.; Cerelli M.J.; Daniloff G.Y.; Kung V.T.  
 CORPORATE SOURCE: Metra Biosystems, 265 N. Wisman Rd., Mountain View, CA 94043-3911, United States.  
 SOURCE: Clinical Chemistry, (1995), 41/11 (1560-1566)  
 CODEN: CLCHAU ISSN: 0009-9147  
 DOCUMENT TYPE: Journal; Article  
 COUNTRY: United States  
 LANGUAGE: English  
 SUMMARY LANGUAGE: English  
 AN 1995:25337878 BIOTECHNO  
 AB Alkaline phosphatase (ALP) is present in human **serum** in the form of several isoenzymes. The two major circulating ALP isoenzymes, bone and liver, are difficult to distinguish because they are the products of a single gene and differ only by posttranslational glycosylation. Quantitative measurement of bone ALP (BAP) activity in **serum** can provide an index for the rate of bone formation. Furthermore, increased BAP activity in **serum** is indicative of bone disorders. We describe a method in which **serum** samples are added to a microtiter plate coated with monoclonal anti-BAP antibody and incubated 3 h at room temperature. After the unbound materials are washed off, the bound BAP activity is measured by adding p-nitrophenyl phosphate substrate. The assay demonstrated no cross-reactivity to intestinal or placental ALP and only 3-8% cross-reactivity to liver ALP. The intraassay (n = 21) CVs were 3.9-5.9%, and interassay (n = 8) CVs were 4.4-7.0%. Comparisons of the assay (y) with an IRMA (x) and a wheat germ agglutinin precipitation method (x') gave regression equations of  $y = 1.32x - 6.4$ ,  $r = 0.99$ , and  $y = 1.41x' + 4.8$ ,  $r = 0.99$ . The assay detected increased BAP in **sera** from patients with osteoporosis, Paget disease, osteomalacia, or primary hyperparathyroidism.

L96 ANSWER 11 OF 32 BIOTECHNO COPYRIGHT 2002 Elsevier Science B.V.DUPLICATE  
 ACCESSION NUMBER: 1995:25068653 BIOTECHNO  
 TITLE: **Spectrophotometric** determination of **alkaline phosphatase** and **.alpha.-fetoprotein** in human **serum** with 5,10,15,20-tetrakis (4-phosphonooxyphenyl)porphine  
 AUTHOR: Kawakami T.; Igarashi S.  
 CORPORATE SOURCE: Department of Materials Science, Faculty of Engineering, Ibaraki University, Nakanarusawa, Hitachi 316, Japan.  
 SOURCE: Analyst, (1995), 120/2 (539-542)  
 CODEN: ANALAO ISSN: 0003-2654  
 DOCUMENT TYPE: Journal; Article  
 COUNTRY: United Kingdom  
 LANGUAGE: English  
 SUMMARY LANGUAGE: English

AN 1995:25068653 BIOTECHNO  
 AB A highly sensitive kinetic **spectrophotometric** determination for alkaline phosphatase (ALP) using 5,10,15,20-tetrakis (4-phosphonooxyphenyl)porphine (approximately  $\Delta \epsilon \cdot \text{simeq} \cdot 3.0 \times 10^4 \text{ l} \cdot \text{mol}^{-1} \cdot \text{cm}^{-1}$ ) was studied. With a measuring time of 30 min, the determination range was  $2.0 \times 10^{-8} \text{ mol l}^{-1}$  to  $8.0 \times 10^{-7} \text{ mol l}^{-1}$  and the detection limit (3s) was  $1.0 \times 10^{-8} \text{ mol l}^{-1}$ . The detection limit became lower on increasing the measuring time. For example, the detection limit was up to  $10^{-7} \text{ mol l}^{-1}$  for 12 h. This method was applied to the assay of ALP in human **serum**. Using 0.50 ml of sample without any pre-treatment, a good correlation was obtained with the p-nitrophenyl phosphate method, and the relative standard deviation

was 4.7%. The method was also applied to the enzyme immunoassay of .alpha.-fetoprotein in human ~~serum~~ and provided a good correlation with the 4-methylumbelliferone derivative ( $r = 0.999$ ) with approximately half the relative standard deviations (1.3-2.7%) obtained with a commercial kit. The proposed method requires very small sample volumes, is highly sensitive, highly accurate, and shows good stability of the reagent; it will be applied to various assays using ALP as a label.

L96 ANSWER 12 OF 32 MEDLINE  
 ACCESSION NUMBER: 92174307 MEDLINE  
 DOCUMENT NUMBER: 92174307 PubMed ID: 1541007  
 TITLE: Simultaneous determinations of liver- and bone-type alkaline phosphatase by curve-fitting of inhibition kinetic data. I. Development and evaluation of an absorbance-based method.  
 AUTHOR: Fitzpatrick C P; Pardue H L  
 CORPORATE SOURCE: Department of Chemistry, Purdue University, West Lafayette, IN 47907-1393.  
 CONTRACT NUMBER: GMS 13326-22,23 (NIGMS)  
 SOURCE: CLINICAL CHEMISTRY, (1992 Feb) 38 (2) 238-46.  
 Journal code: 9421549. ISSN: 0009-9147.  
 PUB. COUNTRY: United States  
 DOCUMENT TYPE: Journal; Article; (JOURNAL ARTICLE)  
 LANGUAGE: English  
 FILE SEGMENT: Priority Journals  
 ENTRY MONTH: 199204  
 ENTRY DATE: Entered STN: 19920424  
 Last Updated on STN: 19970203  
 Entered Medline: 19920403

AB We describe an approach for the simultaneous determination of isoenzymes of alkaline phosphatase (EC 3.1.3.1) based on the kinetic behavior of inhibition reactions. Data for absorbance vs time, collected while enzymes are being inhibited, are fitted with suitable models to obtain results related to activities of the individual isoenzymes. The primary focus is on two-component mixtures of the bone and liver isoenzymes of alkaline phosphatase, but some results are reported for three- and four-component mixtures. Factors studied include choices of inhibitors, buffers, pH, ionic strength, substrate concentration, kinetic models, data ranges, data densities, and data-processing approaches and programs. Criteria used to select optimal conditions include measurement times, detection limits, useful range, and agreement between expected and computed results for mixtures of isoenzymes. For two-component mixtures, a linear least-squares fit of isoenzyme content computed with the curve-fitting method (y) v a comparison method (x) gave  $y = 0.96 (+/- 0.05)x + 3.8 (+/- 3)\%$  with  $r = 0.97$  and standard error of the estimate of 9.4% for a range from 15 to 300 U/L. The pooled relative standard deviation (CV) for results was about 5%. Results were degraded for three- and four-component samples.

L96 ANSWER 13 OF 32 WPIDS (C) 2002 THOMSON DERWENT  
 ACCESSION NUMBER: 1991-143529 [20] WPIDS  
 DOC. NO. NON-CPI: N1991-110372  
 DOC. NO. CPI: C1991-061740  
 TITLE: Method of discriminating human blood and animal blood - using monoclonal antibody specific to human **haemoglobin**.  
 DERWENT CLASS: B04 S03  
 PATENT ASSIGNEE(S): (ASAK) ASAHI BREWERIES LTD  
 COUNTRY COUNT: 1  
 PATENT INFORMATION:

PATENT NO	KIND	DATE	WEEK	LA	PG
JP 03078655	A	19910403	(199120)*		

## APPLICATION DETAILS:

PATENT NO	KIND	APPLICATION	DATE
JP 03078655	A	JP 1989-214922	19890823

PRIORITY APPLN. INFO: JP 1989-214922 19890823

AN 1991-143529 [20] WPIDS

AB JP 03078655 A UPAB: 19930928

A method of discriminating human blood and animal blood comprises judging the origin of sample blood about blood, blood-stain or sample contg. blood by using monoclonal antibody specific to human **hemoglobin**. A reagent for discriminating human blood and animal blood comprising using monoclonal antibody to human **hemoglobin** is also claimed. The monoclonal antibody specific to human **hemoglobin** which is obtd. as follows. Antibody forming cell obtd. by immunising other animals than human being with human **hemoglobin** (antigen) is fused with myeloma cell, resulting hybridoma is cultured in a medium or transplanted into the abdominal cavity of mouse, and screening and cloning are repeated as desired, obtaining anti-human **hemoglobin** monoclonal antibody capable of binding selectively with human **hemoglobin**.

The amt. of sample blood used for the discrimination can be extremely little. The discrimination can be carried out by e.g. ELISA method in which the monoclonal antibody is bound with a sample **blood**, treated with **alkaline phosphatase**-labelled anti-mouse IgG, treated with **p-nitrophenyl phosphate**, and the absorbence of the reaction liq. is measured. By the method, human blood and blood of Japanese monkey, chimpanzee, gorilla, etc. can be discriminated.

USE/ADVANTAGE - The method and the reagent are useful for discriminating the origin of sample blood, human being or other animals, about a sample contg. blood. The discrimination of the origin of blood can be simply and accurately carried out without complex operation such as absorbing operation by utilising high specificity of the monoclonal antibody to human **hemoglobin**. @ (6pp Dwg.No 0/0)@

L96 ANSWER 14 OF 32 WPIDS (C) 2002 THOMSON DERWENT

ACCESSION NUMBER: 1989-218793 [30] WPIDS

DOC. NO. NON-CPI: N1989-166548

DOC. NO. CPI: C1989-097501

TITLE: Monitoring of elimination of gypsum on equipment surfaces - by **spectrophotometric** methods involving determin. of nitro-tri methyl phosphonic acid from change in hydrolysis rate of **nitrophenyl phosphate**.

DERWENT CLASS: E11 J04 M25 S03

INVENTOR(S): DYATLOVA, N M; TEREKHIN, S N; VOLKOVA, N A

PATENT ASSIGNEE(S): (CHRE-R) CHEM REAGENTS PRS; (MOYU) MOSCOW LOMONOSOV UNIV;  
(USTK-R) UST-KAMENOGORSK LEAD

COUNTRY COUNT: 1

PATENT INFORMATION:

PATENT NO	KIND	DATE	WEEK	LA	PG

SU 1434366 A 19881030 (198930)\* 2

## APPLICATION DETAILS:

PATENT NO	KIND	APPLICATION	DATE
SU 1434366	A	SU 1986-4121060	19860707

PRIORITY APPLN. INFO: SU 1986-4121060 19860707

AN 1989-218793 [30] WPIDS

AB SU 1434366 A UPAB: 19930923

The process involves determin. of free nitrilo-trimethyl-phosphonic acid by **spectrophotometric** means, according to the variation in the rate of reaction of the hydrolysis of para-**nitrophenyl phosphate** by **alkaline phosphatase** at pH 10-11.

The method enables the extent to which the gypsum deposits are eliminated to be quantitatively determined.

USE - In non-ferrous metals extn. and refining industries, for removing gypsum deposits from surfaces, partic. in the prodn. of zinc. Bul.40/30.10.88

0/0

L96 ANSWER 15 OF 32 WPIDS (C) 2002 THOMSON DERWENT

ACCESSION NUMBER: 1988-058205 [09] WPIDS

DOC. NO. NON-CPI: N1988-044224

DOC. NO. CPI: C1988-025942

TITLE: Dry analytical element for theophylline determin. - contg. increased concn. of **alkaline phosphatase** isoenzyme to reduce interferences.

DERWENT CLASS: B04 D16 S03

INVENTOR(S): FRICKEY, P H; NORTON, G E

PATENT ASSIGNEE(S): (EAST) EASTMAN KODAK CO

COUNTRY COUNT: 9

PATENT INFORMATION:

PATENT NO	KIND	DATE	WEEK	LA	PG
EP 258035	A	19880302	(198809)*	EN	9
R: CH DE FR GB LI					
JP 63079600	A	19880409	(198820)		
US 4806470	A	19890221	(198910)		7
CA 1290662	C	19911015	(199150)		
EP 258035	B1	19920603	(199223)	EN	10
R: CH DE FR GB LI					
DE 3779519	G	19920709	(199229)		
JP 06098033	B2	19941207	(199502)		6

## APPLICATION DETAILS:

PATENT NO	KIND	APPLICATION	DATE
EP 258035	A	EP 1987-307498	19870825
JP 63079600	A	JP 1987-208367	19870824
US 4806470	A	US 1986-900068	19860825
EP 258035	B1	EP 1987-307498	19870825
DE 3779519	G	DE 1987-3779519	19870803
		EP 1987-307498	19870803
JP 06098033	B2	JP 1987-208367	19870824

## FILING DETAILS:

PATENT NO	KIND	PATENT NO
DE 3779519	G Based on	EP 258035
JP 06098033	B2 Based on	JP 63079600

PRIORITY APPLN. INFO: US 1986-900068 19860825

AN 1988-058205 [09] WPIDS

AB EP 258035 A UPAB: 19930923

Element comprises an absorbent carrier material, a buffer which maintains the pH at 9 or less during the detn. and, in fluid contact, first and second zones, the first zone contg. an isoenzyme of **alkaline phosphatase** which is capable of acting on a substrate for the isoenzyme at a pH of 9 or less and the second zone contg. a substrate for the isoenzyme. The isoenzyme is present in an amt. of at least 100IU/m2.

Pref. the substrate is an organic mono- or diester of phosphoric acid, esp. **p-nitrophenyl phosphate** or 4-(4-nitro-2-methylsulphonyl phenylazo)naphthol-1-phosphate. Pref. the isoenzyme is bovine liver **alkaline phosphatase**.

USE/ADVANTAGE - The element provides a simple and rapid assay for theophylline esp. in body fluids, having all of the advantages of the element in EP184,437. In addn., by increasing the amt. of **alkaline phosphatase** isoenzyme in the element, the assay exhibits less interferences from endogenous phosphatase and **haemoglobin** in the biological test fluids.

0/0

ABEQ DE 3779519 G UPAB: 19930923

Element comprises an absorbent carrier material, a buffer which maintains the pH at 9 or less during the detn. and, in fluid contact, first and second zones, the first zone contg. an isoenzyme of **alkaline phosphatase** which is capable of acting on a substrate for the isoenzyme at a pH of 9 or less and the second zone contg. a substrate for the isoenzyme. The isoenzyme is present in an amt. of at least 100IU/m2.

Pref. the substrate is an organic mono- or diester of phosphoric acid, esp. **p-nitrophenyl phosphate** or 4-(4-nitro-2-methylsulphonyl phenylazo)naphthol-1-phosphate. Pref. the isoenzyme is bovine liver **alkaline phosphatase**.

USE/ADVANTAGE - The element provides a simple and rapid assay for theophylline esp. in body fluids, having all of the advantages of the element in EP184,437. In addn., by increasing the amt. of **alkaline phosphatase** isoenzyme in the element, the assay exhibits less interferences from endogenous phosphatase and **haemoglobin** in the biological test fluids.

ABEQ EP 258035 B UPAB: 19930923

A dry analytical element for the determination of theophylline comprising an absorbent carrier material, a buffer which maintains the pH at 9 or less during the determination and, in fluid contact, first and second zones, the first zone containing an isoenzyme of **alkaline phosphatase** which is capable of acting on a substrate for the isoenzyme at a pH of 9 or less, and the second zone containing a substrate for the isoenzyme, the element characterized in that the **alkaline phosphatase** isoenzyme is present in an amount of at least about 100 I.U./m2.

ABEQ US 4806470 A UPAB: 19930923

Analytical test strip for theophylline comprises a dry absorbent carrier impregnated with a buffer agent (to maintain pH not above 9.0); having two distinct zones, one impregnated with **alkaline phosphatase** isoenzyme (at least 100 int. units/m2) and the other impregnated with a substrate for the isoenzyme, e.g. 4-

**nitrophenyl phosphate** of 4-(2-methylsulphonyl-4-nitrophenyl) azonaphthyl-1-phosphate.

USE - In aq. media, the isoenzyme and organic phosphate react to give a highly coloured prod., but the enzyme is inhibited when theophylline is present in the aq. medium, with no colour formation.

L96 ANSWER 16 OF 32 MEDLINE

ACCESSION NUMBER: 88022450 MEDLINE

DOCUMENT NUMBER: 88022450 PubMed ID: 3662403

TITLE: The effect of haemolysis on the measurement of plasma alkaline phosphatase activity.

AUTHOR: Grosset A; Knapp M L; Mayne P D

CORPORATE SOURCE: Department of Chemical Pathology, Charing Cross and Westminster Medical School, Westminster Hospital, London, UK.

SOURCE: ANNALS OF CLINICAL BIOCHEMISTRY, (1987 Sep) 24 ( Pt 5) 513-7.

Journal code: 0324055. ISSN: 0004-5632.

PUB. COUNTRY: ENGLAND: United Kingdom

DOCUMENT TYPE: Journal; Article; (JOURNAL ARTICLE)

LANGUAGE: English

FILE SEGMENT: Priority Journals

ENTRY MONTH: 198711

ENTRY DATE: Entered STN: 19900305

Last Updated on STN: 19900305

Entered Medline: 19871120

AB The non-dialysable fraction of haemolysate causes an apparent reduction of plasma alkaline phosphatase (ALP) activity using 4-nitrophenylphosphate as substrate. Analyses using four different buffers showed that the decrease in enzyme activity is affected by the buffer used. The percentage reduction in ALP activity is dependent on the initial ALP activity but not on the isoenzyme present. When diethanolamine was used as buffer, sample blanking almost completely compensated for the apparent reduction in enzyme activity. However, when aminomethylpropanol, aminomethylpropanediol and tris-carbonate buffers were used, it appeared that haemolysate reduced the catalytic activity of the enzyme, since sample blank correction had minimal effect on the results.

L96 ANSWER 17 OF 32 CAPLUS COPYRIGHT 2002 ACS

ACCESSION NUMBER: 1986:221645 CAPLUS

DOCUMENT NUMBER: 104:221645

TITLE: Stable alkaline phosphatase color test elements

INVENTOR(S): Nagashima, Minoru; Azuma, Masayuki; Noguchi, Sadao

PATENT ASSIGNEE(S): Kyowa Hakko Kogyo Co., Ltd., Japan

SOURCE: Jpn. Kokai Tokkyo Koho, 5 pp.

CODEN: JKXXAF

DOCUMENT TYPE: Patent

LANGUAGE: Japanese

FAMILY ACC. NUM. COUNT: 1

PATENT INFORMATION:

PATENT NO.	KIND	DATE	APPLICATION NO.	DATE
JP 61015690	A2	19860123	JP 1984-134534	19840629

AB A multilayer test element is developed for alk. phosphatase (I) detn. in clin. specimens. The element consists of (1) a support, (2) a development layer contg. dry microparticles of a I substrate that develops a color upon contact with I, (3) a hydrophilic reagent layer, and (4) optionally a protective film layer. Because of the use of dry microparticles of



substrate, the product is highly stable and can be kept at 5.degree. for >3 mo. Thus, one side of a polyethylene terephthalate film was coated successively with reagent 1 contg. gelatin, Na triisopropyl naphthalenesulfonate, and NaHCO<sub>3</sub>, reagent 2 contg. gelatin, Na triisooctyl naphthalenesulfonate, and 1,2-bis(vinylsulfonyl)ethane, a development layer compn. contg. powd. filter paper, Triton X-100, styrene-glycidyl methacrylate copolymer, MgSO<sub>4</sub>.7H<sub>2</sub>O, and p-nitrophenyl phosphate. The element was contacted with human serum and its optical d. was detd. at 410 nm.

L96 ANSWER 18 OF 32 WPIDS (C) 2002 THOMSON DERWENT  
 ACCESSION NUMBER: 1986-191159 [30] WPIDS  
 DOC. NO. NON-CPI: N1986-142851  
 DOC. NO. CPI: C1986-082319  
 TITLE: Analytical element for theophylline determn. - by determining inhibition of alkaline phosphatase by theophylline at low pH.  
 DERWENT CLASS: B04 D16 J04 S03  
 INVENTOR(S): NORTON, G E  
 PATENT ASSIGNEE(S): (EAST) EASTMAN KODAK CO  
 COUNTRY COUNT: 10  
 PATENT INFORMATION:

PATENT NO	KIND	DATE	WEEK	LA	PG
EP 188372	A	19860723	(198630)*	EN	25
R: CH DE FR GB IT LI					
JP 61170400	A	19860801	(198637)		
US 4782016	A	19881101	(198846)		
CA 1251385	A	19890321	(198916)		
EP 188372	B	19920102	(199202)		
R: CH DE FR GB IT LI					
DE 3683165	G	19920213	(199208)		
JP 06087799	B2	19941109	(199443)		8

## APPLICATION DETAILS:

PATENT NO	KIND	APPLICATION	DATE
EP 188372	A	EP 1986-300226	19860115
JP 61170400	A	JP 1986-6587	19860117
US 4782016	A	US 1985-692473	19850118
JP 06087799	B2	JP 1986-6587	19860117

## FILING DETAILS:

PATENT NO	KIND	PATENT NO
JP 06087799	B2 Based on	JP 61170400

PRIORITY APPLN. INFO: US 1985-692473 19850118

AN 1986-191159 [30] WPIDS

AB EP 188372 A UPAB: 19930922

A dry analytical element for the determn. of theophylline in a human biological fluid comprises an absorbent carrier material (I), a buffer which maintains the pH at 9 or less during the determination and, in fluid contact, first and second zones. The first zone contains an isoenzyme of alkaline phosphatase which is capable of acting on a substrate for the isoenzyme at a pH of 9 or less and the second zone contains a substrate for the isoenzyme. Suitable substrates are organic mono- or diesters of

phosphoric acid e.g. **p-nitrophenyl phosphate** or 4-(4-nitro -2-methyl-sulphonyl phenylazo)naphthol -1-phosphate. The method can also be used in soln..

**ADVANTAGE** - The element provides a rapid and simple assay which avoids the effect of endogenous **alkaline phosphatase** in **serum** samples. The **alkaline phosphatase** does not have significant activity at pH 9 or less.  
0/0

ABEQ EP 188372 B UPAB: 19930922

A dry analytical element for the determination of theophylline in a human biological fluid comprising an absorbent carrier material, a buffer which maintains the pH at 9 or less during the determination and, first and second zones which are in fluid contact, the element characterised in that the first zone contains a nonhuman isoenzyme of alkaline phosphatase which is capable of acting on a substrate for the isoenzyme at a pH of 9 or less, and the second zone contains a substrate for the isoenzyme.

ABEQ US 4782016 A UPAB: 19930922

Dry analytical element for theophylline determ. in human biological fluid comprises an absorbent carrier, buffer to maintain pH 9 or less and, in fluid contact, a first zone contg. isoenzyme of alkaline phosphatase which can act on a substrate for the isoenzyme at pH 9 or less and a second zone contg. substrate for the isoenzyme. The substrate may be organo mono- or di-ester of phosphoric acid and the isoenzyme may be bovine liver alkaline phosphatase.

**ADVANTAGE** - Simple, rapid analysis can be effected without laborious pretreatment or extn. techniques.

L96 ANSWER 19 OF 32 BIOTECHNO COPYRIGHT 2002 Elsevier Science B.V.

ACCESSION NUMBER: 1986:16050862 BIOTECHNO

TITLE: Elevation of **alkaline phosphatase** by retinol in bovine endothelial cells and its possible relationship to lipid biosynthesis

AUTHOR: Adams S.E.; Bishop E.J.; Melnykovych G.

CORPORATE SOURCE: Department of Microbiology, University of Kansas Medical Center, Kansas City, KS 66103, United States.

SOURCE: Biochimica et Biophysica Acta - Molecular Cell Research, (1986), 885/2 (146-153)

CODEN: BAMRDP

DOCUMENT TYPE: Journal; Article

COUNTRY: Netherlands

LANGUAGE: English

AN 1986:16050862 BIOTECHNO

AB Alkaline phosphatase (EC 3.1.3.1) activity in bovine aortic endothelial cells in culture was stimulated in a synergistic manner by 10<sup>sup.6</sup> M retinol and by 10<sup>sup.7</sup> M dexamethasone. An early exposure to retinol was required for maximum stimulation and could be reproduced by the addition, during growth, of 2  $\mu$ g/ml compactin. The induced enzyme activity in cell lysates prepared from cells treated with retinol and dexamethasone had a V(max) that was 50-fold that of the controls. The stimulatory effect of retinol could be partially reversed by the addition of sonic dispersions made from cholesterol and phosphatidylcholine. The incorporation of <sup>14</sup>C-acetate into saponifiable and non-saponifiable cellular lipids was inhibited by 10<sup>sup.6</sup> M retinol but the activities of 3-hydroxy-3-methylglutaryl coenzyme A reductase (EC 1.1.1.34) and 3-hydroxy-3-methylglutaryl coenzyme A synthase (EC 4.1.3.5) remained unaffected. The results suggest that retinol might inhibit lipid biosynthesis through an alternate mechanism.

L96 ANSWER 20 OF 32 MEDLINE

ACCESSION NUMBER: 83180697 MEDLINE

DOCUMENT NUMBER: 83180697 PubMed ID: 6404566  
TITLE: A reference method for measurement of alkaline phosphatase activity in human serum.  
AUTHOR: Tietz N W; Burtis C A; Duncan P; Ervin K; Petittclerc C J; Rinker A D; Shuey D; Zygowicz E R  
SOURCE: CLINICAL CHEMISTRY, (1983 May) 29 (5) 751-61.  
Journal code: 9421549. ISSN: 0009-9147.  
PUB. COUNTRY: United States  
DOCUMENT TYPE: Journal; Article; (JOURNAL ARTICLE)  
LANGUAGE: English  
FILE SEGMENT: Priority Journals  
ENTRY MONTH: 198306  
ENTRY DATE: Entered STN: 19900318  
Last Updated on STN: 19900318  
Entered Medline: 19830610

AB We present an official AACC reference method for the measurement of alkaline phosphatase, the culmination of optimization experiments conducted by a group of independent laboratories. The details of this method and evaluation of factors affecting the measurement are described. A metal ion buffer has been incorporated that maintains optimal and constant concentrations of zinc(II) and magnesium(II) ions. Final reaction conditions are: pH (30 degrees C), 10.40 +/- 0.05; 2-amino-2-methyl-1-propanol buffer, 0.35 mol/L; **4-nitrophenyl phosphate**, 16.0 mmol/L; magnesium acetate, 2.0 mmol/L; zinc sulfate, 1.0 mmol/L; and N-(2-hydroxyethyl)ethylenediaminetriacetic acid, 2.0 mmol/L.

L96 ANSWER 21 OF 32 BIOSIS COPYRIGHT 2002 BIOLOGICAL ABSTRACTS INC.DUPLICATE  
4

ACCESSION NUMBER: 1983:70660 BIOSIS  
DOCUMENT NUMBER: BR24:70660  
TITLE: L PHENYL ALANINE INHIBITION OF HUMAN **ALKALINE PHOSPHATASES** EC-3.1.3.1 WITH P **NITROPHENYL PHOSPHATE** AS SUBSTRATE.  
AUTHOR(S): KOMODA T; HOKARI S; SONODA M; SAKAGISHI Y; TAMURA T  
CORPORATE SOURCE: DEP. BIOCHEM., SAITAMA MED. SCH., 38 MOROYAMA, IRUMA-GUN, SAITAMA 350-04, JPN.  
SOURCE: Clin. Chem. (Winston-Salem, N. C.), (1982) 28 (12), 2426-2428.  
CODEN: CLCHAU. ISSN: 0009-9147.  
FILE SEGMENT: BR; OLD  
LANGUAGE: English

L96 ANSWER 22 OF 32 EMBASE COPYRIGHT 2002 ELSEVIER SCI. B.V.  
ACCESSION NUMBER: 82127132 EMBASE  
DOCUMENT NUMBER: 1982127132  
TITLE: Reactivity and stability of microencapsulated placental **alkaline phosphatase**.  
AUTHOR: Takenaka H.; Kawashima Y.; Chikamatsu Y.; Ando Y.  
CORPORATE SOURCE: Gifu Coll. Pharm., Mitahora, Gifu, 502, Japan  
SOURCE: Chemical and Pharmaceutical Bulletin, (1982) 30/2 (695-701).  
CODEN: CPBTAL  
COUNTRY: Japan  
DOCUMENT TYPE: Journal  
FILE SEGMENT: 037 Drug Literature Index  
LANGUAGE: English

L96 ANSWER 23 OF 32 BIOSIS COPYRIGHT 2002 BIOLOGICAL ABSTRACTS INC.  
ACCESSION NUMBER: 1983:194316 BIOSIS

DOCUMENT NUMBER: BA75:44316  
TITLE: KINETIC STUDIES OF THE TRANS PHOSPHORYLATION REACTIONS  
CATALYZED BY **ALKALINE PHOSPHATASE** FROM  
ESCHERICHIA-COLI HYDROLYSIS OF P **NITROPHENYL**  
**PHOSPHATE** AND O CARBOXYPHENYL PHOSPHATE IN PRESENCE  
OF TRIS HYDROXYMETHYLAMINO METHANE.  
AUTHOR(S): ROIG M G; BURGUILLO F J; DEL ARCO A; USERO J L; IZQUIERDO  
C; HERRAEZ M A  
CORPORATE SOURCE: DEP. QUIM. FISICA, FAC. QUIM., UNIV. SALAMANCA, SALAMANCA,  
SPAIN.  
SOURCE: INT J BIOCHEM, (1982) 14 (7), 655-666.  
CODEN: IJBOBV. ISSN: 0020-711X.  
FILE SEGMENT: BA; OLD  
LANGUAGE: English  
AB Transphosphorylation of p-nitrophenyl phosphate and o-carboxyphenyl  
phosphate to Tris by E. coli alkaline phosphatase was studied at alkaline  
and acid pH. The rate of release for all reaction products was  
Tris-dependent for both substrates, with a slight maximum for phenol at  
alkaline pH. These dependencies were analyzed from a mechanistic  
standpoint. Individual constants of rate of a simple transphosphorylation  
mechanism were determined. At high Tris concentrations (> 1.0 M) a slight  
competitive inhibition was observed. Inhibition in NH<sub>4</sub><sup>+</sup>-NH<sub>3</sub>Cl buffer was  
found at alkaline pH but not at acid pH. The nonprotonated NH<sub>2</sub> group of  
Tris is probably responsible for inhibition. Complexes are apparently  
formed between Tris and the enzyme. Other possible alternatives are also  
analyzed.

L96 ANSWER 24 OF 32 MEDLINE  
ACCESSION NUMBER: 81210644 MEDLINE  
DOCUMENT NUMBER: 81210644 PubMed ID: 7237780  
TITLE: Response-surface-optimized, zinc-enhanced assay for serum  
alkaline phosphatase.  
AUTHOR: Thompson J C; Hodges C T; Dobler G L; Williamson J A Jr  
SOURCE: CLINICAL CHEMISTRY, (1981 Jul) 27 (7) 1171-5.  
Journal code: 9421549. ISSN: 0009-9147.  
PUB. COUNTRY: United States  
DOCUMENT TYPE: Journal; Article; (JOURNAL ARTICLE)  
LANGUAGE: English  
FILE SEGMENT: Priority Journals  
ENTRY MONTH: 198108  
ENTRY DATE: Entered STN: 19900316  
Last Updated on STN: 19970203  
Entered Medline: 19810827  
AB An optimized assay for alkaline phosphatase (EC 3.1.3.1) is reported. A  
discrete analyzer, the DuPont Automatic Clinical Analyzer (aca), was used  
for this study. The assay is based on results of response-surface  
experimental co-optimization techniques, and response is enhanced over the  
present aca assay. A key feature is the incorporation of zinc ions, both  
to fully optimize the assay and to reduce the sensitivity of measured  
activity to zinc-binding impurities in the buffer, 2-amino-2-methyl-1-  
propanol. In addition, a simple technique is described for measuring  
relative concentrations of zinc-binding impurities in this buffer. These  
features should be considered in the design of any assay for alkaline  
phosphatase that is based on **p-nitrophenyl**  
**phosphate** as substrate and 2-amino-2-methyl-1-propanol as buffer.

L96 ANSWER 25 OF 32 MEDLINE DUPLICATE 5  
ACCESSION NUMBER: 81088720 MEDLINE  
DOCUMENT NUMBER: 81088720 PubMed ID: 7449094  
TITLE: **4-nitrophenyl phosphate**

--characterization of high-purity materials for measuring alkaline phosphatase activity in human serum.

AUTHOR: Bowers G N Jr; McComb R B; Upreti A  
SOURCE: CLINICAL CHEMISTRY, (1981 Jan) 27 (1) 135-43.  
Journal code: 9421549. ISSN: 0009-9147.

PUB. COUNTRY: United States  
DOCUMENT TYPE: Journal; Article; (JOURNAL ARTICLE)  
LANGUAGE: English  
FILE SEGMENT: Priority Journals  
ENTRY MONTH: 198103  
ENTRY DATE: Entered STN: 19900316  
Last Updated on STN: 19900316  
Entered Medline: 19810324

AB We studied 53 lots of **4-nitrophenyl phosphate** (I), obtained from 20 different commercial suppliers, and used this information to set specifications for it. Using these well-defined specifications, we classified 21 lots of I as "unacceptable," 26 lots as "borderline," and six as "acceptable." All lots were shown to contain some 4-nitrophenol and inorganic phosphate. However, "acceptable" I had < 0.3 mmol of 4-nitrophenol and < 10 mmol of inorganic phosphate per mole of I. The mole concentration of I (based on disodium hexahydrate, formula weight 371) was determined by enzymic conversion to 4-nitrophenol in five lots of "acceptable" materials. The mole fraction of I ranged from 0.982 to 0.998. From these measurements and from estimates of impurities that absorb at 311 nm, as determined by liquid chromatography and spectrophotometry at other wavelengths, our best estimate of the molar absorptivity of I at 311 nm in 10 mmol/L NaOH at 25 degrees C was 9867 L x mol<sup>-1</sup> x cm<sup>-1</sup>, with a total uncertainty of 76 L x mol<sup>-1</sup> x cm<sup>-1</sup>. We recommend that I used in clinical laboratories for measurement of alkaline phosphatase activity in serum meet the specifications given in this paper: I content > 98%, maximum activity > 98% in comparative testing with other "acceptable" lots of I, and impurities not to exceed the values cited above.

L96 ANSWER 26 OF 32 BIOTECHNO COPYRIGHT 2002 Elsevier Science B.V.

ACCESSION NUMBER: 1981:11188992 BIOTECHNO

TITLE: **4-Nitrophenyl phosphate**

-characterization of high-purity materials for measuring **alkaline phosphatase** activity in human serum

AUTHOR: Bowers Jr. G.N.; McComb R.B.; Upreti A.

CORPORATE SOURCE: Clin. Chem. Lab., Dept. Pathol., Hartford Hosp., Hartford, Conn. 06115, United States.

SOURCE: Clinical Chemistry, (1981), 27/1 (135-143)

CODEN: CLCHAU

DOCUMENT TYPE: Journal; Article

COUNTRY: United States

LANGUAGE: English

AN 1981:11188992 BIOTECHNO

AB We studied 53 lots of 4-nitrophenyl phosphate (I), obtained from 20 different commercial suppliers, and used this information to set specifications for it. Using these well-defined specifications, we classified 21 lots of I as 'unacceptable,' 26 lots as 'borderline,' and six as 'acceptable.' All lots were shown to contain some 4-nitrophenol and inorganic phosphate. However, 'acceptable' I had <0.3 mmol of 4-nitrophenol and <10 mmol of inorganic phosphate per mole of I. The mole concentration of I (based on disodium hexahydrate, formula weight 371) was determined by enzymic conversion to 4-nitrophenol in five lots of 'acceptable' materials. The mole fraction of I ranged from 0.982 to 0.998. From these measurements and from estimates of impurities that absorb at 311 nm, as determined by liquid chromatography and

spectrophotometry at other wavelengths, our best estimate of the molar absorptivity of I at 311 nm in 10 mmol/L NaOH at 25.degree.C was 9867 L.mol.sup.-.sup.1.cm.sup.-.sup.1, with a total uncertainty of 76 L.mol.sup.-.sup.1.cm.sup.-.sup.1. We recommend that I used in clinical laboratories for measurement of alkaline phosphatase activity in serum meet the specifications given in this paper: I content >98%, maximum activity >98% in comparative testing with other 'acceptable' lots of I, and impurities not to exceed the values cited above.

L96 ANSWER 27 OF 32 MEDLINE  
ACCESSION NUMBER: 80104108 MEDLINE  
DOCUMENT NUMBER: 80104108 PubMed ID: 118558  
TITLE: An enzyme linked immunosorbent assay (ELISA) for detecting IgG sensitized erythrocytes.  
AUTHOR: Bruner K W Jr; Kissling C W  
SOURCE: TRANSFUSION, (1979 Nov-Dec) 19 (6) 773-7.  
Journal code: 0417360. ISSN: 0041-1132.  
PUB. COUNTRY: United States  
DOCUMENT TYPE: Journal; Article; (JOURNAL ARTICLE)  
LANGUAGE: English  
FILE SEGMENT: Priority Journals  
ENTRY MONTH: 198003  
ENTRY DATE: Entered STN: 19900315  
Last Updated on STN: 19900315  
Entered Medline: 19800317

AB This paper describes an Enzyme Linked Immunosorbent Assay (ELISA) for detecting IgG sensitized erythrocytes utilizing a commercially available anti-human IgG conjugated with alkaline phosphatase. Erythrocyte hemolysis in the assay was minimized by dissolving the p-nitrophenyl phosphate substrate in a carbonate-bicarbonate buffer. Nonspecific absorption of the enzyme conjugate to erythrocytes and glassware was reduced by adding 1% bovine serum albumin to wash solutions. Assay sensitivity was increased with greater concentrations of enzyme conjugate and erythrocytes in the incubation stage. The sensitivity of the described ELISA procedure is approximately equal to that of the standard antiglobulin test. Some possible future applications of ELISA in the blood bank are discussed.

L96 ANSWER 28 OF 32 BIOSIS COPYRIGHT 2002 BIOLOGICAL ABSTRACTS INC.

ACCESSION NUMBER: 1978:195541 BIOSIS  
DOCUMENT NUMBER: BA66:8038  
TITLE: 4 NITRO PHENOL IN 4 NITROPHENYL PHOSPHATE  
A SUBSTRATE FOR ALKALINE PHOSPHATASE AS  
MEASURED BY PAIRED ION HIGH PERFORMANCE LIQUID  
CHROMATOGRAPHY.  
AUTHOR(S): CULBRETH P H; DUNCAN I W; BURTIS C A  
CORPORATE SOURCE: BUR. LAB., CENT. DIS. CONTROL, PUBLIC HEALTH SERV.,  
ATLANTA, GA. 30333, USA.  
SOURCE: CLIN CHEM, (1977) 23 (12), 2288-2291.  
CODEN: CLCHAU. ISSN: 0009-9147.  
FILE SEGMENT: BA; OLD  
LANGUAGE: English

AB Paired-ion high-performance liquid chromatography was used to determine the 4-nitrophenol content of 4-nitrophenyl phosphate, a substrate for alkaline phosphatase analysis. This was done on a reversed-phase column with a mobile phase of methanol/water, 45/55 by vol, containing 3 ml of tetrabutylammonium phosphate reagent per 200 ml of solvent. At a flow rate of 1 ml/min, 4-nitrophenol was eluted at 9 min and monitored at 404 nm; 4-nitrophenyl phosphate was eluted at 5 min and could be monitored at 311 nm. Samples of 4-nitrophenyl phosphate obtained from several sources

contained 0.3-7.8 mol of 4-nitrophenol per mol of 4-nitrophenyl phosphate.

L96 ANSWER 29 OF 32 BIOSIS COPYRIGHT 2002 BIOLOGICAL ABSTRACTS INC.

ACCESSION NUMBER: 1978:144420 BIOSIS

DOCUMENT NUMBER: BA65:31420

TITLE: TEMPERATURE DEPENDENCE OF THE ABSORBANCE OF ALKALINE SOLUTIONS OF 4-NITROPHENYL PHOSPHATE A POTENTIAL SOURCE OF ERROR IN THE MEASUREMENT OF ALKALINE PHOSPHATASE ACTIVITY.

AUTHOR(S): BURTIS C A; SEIBERT L E; BAIRD M A; SAMPSON E J

CORPORATE SOURCE: BUR. LAB., CENT. DIS. CONTROL, ATLANTA, GA. 30333, USA.

SOURCE: CLIN CHEM, (1977) 23 (9), 1541-1547.

CODEN: CLCHAU. ISSN: 0009-9147.

FILE SEGMENT: BA; OLD

LANGUAGE: English

AB The absorbance of an alkaline solution of 4-nitrophenyl phosphate was a function of temperature. Quantitative evaluation of this phenomenon indicated that it depended on the concentration of the compound and was independent of source, buffer concentration and pH above 9.0; was reversible; was not a result of alkaline hydrolysis or 4-nitrophenol contamination; and correlated with a temperature-induced shift of its absorbance spectrum. The phenomenon may represent a potential analytical problem in methods for alkaline phosphatase in which this compound was the substrate. If thermal equilibrium was not reached and maintained during an alkaline phosphatase assay, the thermochromic response would be included in the measured rate. The magnitude of this error depended on the thermal response and control characteristics of each particular instrument and the reaction conditions under which such an analysis was performed.

L96 ANSWER 30 OF 32 MEDLINE

ACCESSION NUMBER: 75218962 MEDLINE

DOCUMENT NUMBER: 75218962 PubMed ID: 1153923

TITLE: A continuous-flow method for the determination of the activity of serum alkaline phosphatase in diethanolamine buffer.

AUTHOR: Viitala A J; Jokela H A; Penttila I M; Nummi S

SOURCE: SCANDINAVIAN JOURNAL OF CLINICAL AND LABORATORY

INVESTIGATION, (1975 May) 35 (3) 267-73.

Journal code: 0404375. ISSN: 0036-5513.

PUB. COUNTRY: Norway

DOCUMENT TYPE: Journal; Article; (JOURNAL ARTICLE)

LANGUAGE: English

FILE SEGMENT: Priority Journals

ENTRY MONTH: 197511

ENTRY DATE: Entered STN: 19900310

Last Updated on STN: 19900310

Entered Medline: 19751108

AB A procedure for determination of serum alkaline phosphatase activity (EC 3.1.3.1) in diethanolamine (DEA) buffer with an AutoAnalyzer II apparatus was designed. The buffer used was 1.0 mol/l DEA-HC buffer, pH 9.8 at 37 degree C, containing 0.5 mmol/l of MgCl<sub>2</sub> and 10 mmol/l of substrate **4-nitrophenyl-phosphate**. The reaction time was about 3 min at 37 degree C. The enzyme activity (U/l) was calculated by determining the amount of 4-nitrophenol formed in reaction. A sampling rate of 70 samples per hour can be used with good linearity up to 1000 U/l. The results obtained by the new continuous-flow system were compared with those measured by the kinetic method according to the Scandinavian recommendation (10). A close correlation between the two methods was observed.

L96 ANSWER 31 OF 32 MEDLINE  
 ACCESSION NUMBER: 75050823 MEDLINE  
 DOCUMENT NUMBER: 75050823 PubMed ID: 4428830  
 TITLE: [Hydrolysis of arylphosphatates by multiple forms of  
 alkaline phosphatase. Studies on human alkaline  
 phosphatases, II (author's transl)].  
 Hydrolyse von Arylphosphaten durch multiple Formen  
 alkalischer Phosphatasen. Untersuchungen uber alkalische  
 Phosphatasen menschlicher Gewebe, II.  
 AUTHOR: Lorentz K; Flatter B; Voss J; Heydrich D  
 SOURCE: ZEITSCHRIFT FUR KLINISCHE CHEMIE UND KLINISCHE BIOCHEMIE,  
 (1974 Feb) 12 (2) 87-91.  
 Journal code: 7607250. ISSN: 0044-2933.  
 PUB. COUNTRY: GERMANY, WEST: Germany, Federal Republic of  
 DOCUMENT TYPE: Journal; Article; (JOURNAL ARTICLE)  
 LANGUAGE: German  
 FILE SEGMENT: Priority Journals  
 ENTRY MONTH: 197502  
 ENTRY DATE: Entered STN: 19900310  
 Last Updated on STN: 19970203  
 Entered Medline: 19750201

L96 ANSWER 32 OF 32 WPIDS (C) 2002 THOMSON DERWENT  
 ACCESSION NUMBER: 1973-28960U [21] WPIDS  
 TITLE: Alkaline phosphatase determination - using 4-  
**nitrophenyl phosphate** as substrate and  
 diethanolamine organic salts as buffer.  
 DERWENT CLASS: B04 S03 S05  
 PATENT ASSIGNEE(S): (HAUR) HAURY CHEMISCHE FAB H  
 COUNTRY COUNT: 1  
 PATENT INFORMATION:

PATENT NO	KIND	DATE	WEEK	LA	PG
DE 2154993	A		(197321)*		

PRIORITY APPLN. INFO: DE 1971-2154993 19711105

AN 1973-28960U [21] WPIDS

AB DE 2154993 A UPAB: 19930831

In the determination of **alkaline phosphatase** in body fluids (esp. **serum**) using **p-nitrophenyl phosphate** (I) as substrate, the p-nitrophenol released by the enzymatic reaction being determined photometrically, an organic salt of diethanolamine (pref. the formate or glutarate) is used instead of diethanolamine hydrochloride as buffer substance. This gives increased turnover of substrate and increases the sensitivity of the determination. Since solutions of (I) are unstable, it is advantageous to prepare a lyophilisate of (I) together with a protective substance which does not interfere with the enzymatic reaction (pref. glucose) and to add the lyophilisate to the buffer soln. immediately before carrying out the determination.